nenor[®]

Sereno

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

Wersja: 1.1


Model: YM-2T9

Table of contents

Spis treści

4	Instrukcja obsługi	PL
16	User manual	EN
27	Bedienungsanleitung	DE
39	Návod k použití	CZ
50	Používateľská príručka	SK
61	Használati utasítás	HU
73	Användarhandbok	SE
84	Käyttäjän käsikirja	FI
95	Brukerhåndbok	NO
106	Brugermanual	DK
118	Gebruikershandleiding	NL
130	Manual del usuario	ES
142	Manuale utente	IT
154	Manual de utilizare	RO
166	Korisnički priručnik	HR/BA
177	Uputstvo za upotrebu	RS/ME
188	Uporabniški priročnik	SI
199	Εγχειρίδιο χρήστη	GR
212	Lietotāja rokasgrāmata	LV
223	Vartotojo vadovas	LT
234	Kasutusjuhend	ET
245	Manuel de l'utilisateur	FR
257	Керівництво користувача	UA
269	Ръководство за употреба	BG
281	Упатство за употреба	MK
294	Manual do utilizador	PT

A**B**

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- a) "Bio-effects and Safety of Diagnostic Ultrasound" issued by AIUM in 1993
- b) "Medical Ultrasound Safety" issued by AIUM in 1994
- c) „Acoustic Output Measurement Standard for Diagnostic Ultrasound Equipment, Revision” issued by AIUM/NEMA in 2004
- d) „Standard for real-time display of thermal and mechanical acoustic output indices on diagnostic ultrasound equipment, Revision 2” issued by AIUM/NEMA in 2004
- e) „Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Probes” issued in 2008.

- f) "Medical electrical equipment—Part2-37: Particular requirements for the basic safety and essential performance of ultrasonic medical diagnostic and monitoring equipment" issued by IEC in 2007.

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USER MANUAL

Terms used in this manual

This manual has been prepared to present key safety rules.

WARNING

The **WARNING** label indicates actions or situations that could result in **injury or death**.

CAUTION

The **CAUTION** label indicates actions or situations that could **damage the device, cause incorrect data or invalidate the procedure**.

NOTE

The **NOTE** label provides useful information about a function or procedure.

Dear Customer,

Thank you for purchasing the Neno Sereno ultrasound foetal heart rate detector. This product is designed for the accurate detection, listening and display of foetal heart rates from early pregnancy to delivery in single pregnancies, both at home and in clinical settings. This user manual contains all the information necessary for proper use. Following the recommendations will maintain the high quality of the device and extend its service life. Before use, please read this manual and keep it for future reference.

01. WARNINGS AND PRECAUTIONS

WARNINGS and PRECAUTIONS must be strictly observed. To avoid the risk of injury and ensure the correct operation of the device, follow the recommendations below:

WARNINGS:

1. The device is not a substitute for a professional examination. If an abnormal foetal heartbeat is detected or cannot be located, the pregnant woman should go to the hospital immediately for a medical consultation. If the woman does not feel the foetus moving, she should also go to a medical facility immediately to seek professional help.
2. The device is not explosion-proof and must not be used in the presence of flammable anaesthetics.
3. Exposure to ultrasound should be kept to a minimum.
4. The device is not protected against defibrillation.
5. Do not use the device simultaneously with high-frequency surgical equipment.
6. Do not use simultaneously with other ultrasound equipment or equipment connected to the patient, such as pacemakers or other electrical devices.
7. Magnetic and electrical fields may interfere with the operation of the device. All nearby devices must comply with the relevant electromagnetic compatibility (EMC) standards. X-ray or magnetic resonance imaging equipment may cause interference.
8. To avoid misdiagnosis: (a) move the probe to obtain the best FHR signal depending on the position of the foetus; (b) if foetal death is suspected, verify the results with another method.
9. If problems occur, stop using the device and contact an authorised service centre.
10. Do not disassemble, repair or modify the device. If you encounter any problems, contact an authorised service centre.
11. Do not heat or throw batteries into fire – this may cause an explosion.

PRECAUTIONS

1. Read the operating instructions and safety precautions carefully before use.
2. Before use, check the equipment for visible signs of damage that could affect staff safety or testing capabilities. If damage is detected, it is recommended that the device be replaced. Always check that the main unit and probe are in good condition before use.
3. Do not touch the patient, power socket and probe at the same time.
4. Only use original parts and accessories included in the kit. Parts and accessories not approved for use with this device may not meet the expected specifications or may damage the device.
5. Do not immerse, throw or soak the detector and battery in water or seawater.

- Electrically powered medical devices must be installed and put into service in accordance with the electromagnetic compatibility information in this user manual.
- Do not connect any devices or accessories to the device that have not been approved by the manufacturer or certified as compliant with IEC 60601-1. The operation or use of unapproved devices or accessories with the device has not been tested or supported, and the operation and safety of the device cannot be guaranteed.
- Do not use the device near other devices or stack it with other devices, and if necessary, observe the device to ensure that it is functioning properly in the configuration in which it will be used.
- Portable and mobile radio communication devices may affect the operation of medical electrical equipment. Avoid strong sources of electromagnetic radiation, such as radio transmitters and mobile phones.
- The device is not MR safe. Do not use this device near or in a magnetic resonance imaging (MRI) environment, as the strong magnetic field may damage the device or pose a safety hazard.
- Do not service or maintain the device or any accessories that are used on a patient.
- Handle the probe with care – do not drop it on hard surfaces.
- Keep the device clean and avoid vibration.
- Do not use strong solvents (e.g. acetone) or abrasives (e.g. steel wool, polishing pastes) for cleaning. Do not leave any liquids on the surface after cleaning.
- Do not immerse the main unit or probe cable in any liquid – the device is not waterproof.
- Do not extend the probe cable beyond 2 metres.
- To avoid interference, do not remove or insert the probe cable while the device is switched on. Remove the probe before switching on and insert it after switching off the device.
- Reduce the volume before using headphones. Ill-fitting headphones may impair sound reception.
- If the patient is sensitive to conductive gel, use water or oil instead of gel.
- The device can only be switched on when the battery cover is closed.
- Do not throw batteries into a fire – they may explode.
- Do not dispose of the batteries or the product in mixed municipal waste. Follow the applicable regulations for the disposal of electronic equipment and batteries.
- If the device will not be used for a long time, remove the batteries to avoid damage.
- Do not insert new and partially used batteries into the device at the same time. This may cause damage.
- Do not use or store the device in areas where it may be exposed to harmful vapours or volatile substances.
- Keep the device out of the reach of infants and children.

02. INTENDED USE AND CONTRAINDICATIONS

The Neno Sereno foetal heart rate detector uses continuous wave (CW) Doppler technology to detect the foetal heartbeat, display its frequency (FHR) and reproduce the sound of the heartbeat during a single pregnancy – from 12 weeks of pregnancy until delivery.

The product is a portable, battery-powered obstetric device intended for use in medical facilities and at home for daily self-monitoring by pregnant women.

No contraindications for the use of the device have been identified. When used correctly, the device does not cause any side effects.

03. EXPLANATION OF SYMBOLS (FIGURE A)

- CE conformity marking with the number of the notified body
- Medical device
- Manufacturer
- Authorised representative in the European Community
- Importer
- Do not dispose of the product in a mixed municipal waste container. Dispose of the product in accordance with the guidelines for the disposal of electronic equipment of this type
- Device with BF parts
- Read the instructions for use
- Date of manufacture
- Unique product identification code
- Production batch number
- Tightness class
- Warning
- Flat cardboard packaging
- Polyethylene terephthalate

16. Storage temperature range
17. Storage humidity range
18. Storage atmospheric pressure range.

04. PRODUCT DESCRIPTION

Device construction – see FIGURE B

1. Main unit
2. OLED display
3. Power/menu button
4. Volume up button
5. Volume decrease button
6. Ultrasonic probe
7. Battery cover
8. Micro USB port for connecting the probe
9. Headphone jack

Display – see FIGURE C

1. Foetal heart rate (FHR). No FHR value, which means that the probe is not connected properly or there is no foetal heart sound
2. Foetal heart sound symbol
3. Speaker volume level
4. Average foetal heart rate (FHR)
5. Measurement time
6. Current volume
7. Battery charge level
8. Frequency of the probe connected to the main unit

05. INSTALLATION

Before use, ensure that the kit is complete and that the individual components show no signs of damage.

Installing the battery

To install the battery, first open the battery compartment cover **FIGURE B7**. Insert the battery, paying attention to the correct polarity, then press the cover and slide it in the opposite direction to removing the battery until you hear a distinctive click. Always ensure that the battery compartment is properly closed before using the detector.

CAUTION: Incorrectly installed batteries may damage the device.

Connecting and operating the probe

Hold the main unit of the detector with one hand and connect the probe connector to the micro USB port with the other hand **FIGURE B8**.

NOTE: Never connect or disconnect the probe while the monitor is turned on – always connect the probe before turning on the device and disconnect it after turning it off.

Using headphones

To connect headphones to the detector, insert the headphone plug into the socket **FIGURE B9**. Once the headphones are connected, the speaker on the device will be automatically muted.

06. INSTRUCTIONS FOR USE (FHR Examination)

To perform a foetal heart rate (FHR) examination correctly, follow the procedure below:

1. Feel for the position of the foetus with your hand to find the best position for detecting the heartbeat.
NOTE: The best signal quality is obtained only when the probe is placed in the optimal position. Avoid areas where loud sounds from the placenta or blood flow in the umbilical cord can be heard. The mother's heart may be detected if the probe is pressed firmly against the mother's blood vessel.
 - If the foetus is in the head-down position and the mother is lying on her back, the heart is usually best heard on **the midline below the umbilicus**.
 - During the examination, avoid lying on your back for long periods of time due to the risk of supine hypotension – **sitting or side** positions are preferred.
 - An accurate FHR reading can only be obtained when an audible foetal heart signal is detected.
2. Turn on the detector using the power button **FIGURE B3**.
3. Apply a small amount of gel to the acoustic surface of the probe.
4. Place the probe on the body at an appropriate angle to obtain the best audio signal. The probe may become slightly warm when placed on the body.
5. Set the appropriate volume using the volume up (**FIGURE B4**) or volume down (**FIGURE B5**) buttons.

- Once the foetal heartbeat signal is detected, you will hear a clear heartbeat sound in the speakers or headphones, and the numerical value of the foetal heart rate (FHR) will appear on the display.
- To switch the measurement result display from numerical values to a graph, briefly press the power button (FIGURE B3). Pressing it again will display the settings screen, and pressing it again will change the view to the numerical value of the measured foetal heart rate.
- After completing the examination, switch off the device, wipe the probe with a clean, soft cloth to remove any gel residue and place it in the holder. There are two ways to switch off the device:
 - Press and hold the **power** button for about 3 seconds (FIGURE B3),
 - The device will turn off automatically after 23 seconds of inactivity.

07. STORAGE, MAINTENANCE, CLEANING AND DISINFECTION

The device is a medical device. Follow your doctor's recommendations and use the device correctly.

STORAGE

The detector should be stored in a dry, clean place, away from sunlight. The device is not a toy – keep it out of the reach of children.

MAINTENANCE

Before each use, check the device and accessories for visible damage and ensure that the device is working properly and is in good working order.

CAUTION: If you detect any damage, **do not use the device**. Contact an authorised service centre.

- Avoid scratches and damage to the display. Sudden changes in temperature or humidity may cause condensation to form on the screen.
- Keep the exterior of the device clean and free of dust and dirt.
- Handle the probe with care to avoid damaging the cover, piezoelectric crystals, or moving mechanism. Do not touch the probe with hard or sharp objects.
- Do not crush or strain the probe cable.
- After use, wipe any gel residue from the probe.

CLEANING

The main unit of the foetal heart rate detector should be wiped with a clean paper towel, soft cloth or damp cloth after each use. If a damp cloth or tissue is used, the monitor should be dried. **CAUTION:** Do not immerse the detector in water or other liquids for cleaning!

The foetal heart rate detector probe should be cleaned after each use as follows:

- Remove gel from the probe using a soft, clean cloth.
- Disconnect the probe from the detector main unit.
- Use a damp soft cloth to remove any residue from the probe or cable. Do not reuse the same cloth.
- If rinsing is necessary, take care not to expose the system connector to moisture or liquids. Mild soap or enzyme-based detergent may be used in accordance with the manufacturer's recommendations.
- After cleaning, do not leave any liquid on the surface. Use a soft, clean, dry, lint-free cloth or tissue to thoroughly dry the probe and cable.

CAUTION: Do not use strong solvents or abrasives for cleaning.

DISINFECTION

Before disinfecting the main unit of the detector, remove the batteries and then clean the screen and housing with a soft cloth soaked in ethyl or isopropyl alcohol (70–90%). During disinfection, avoid excessive wetting of the surface and contact of the liquid with the interior of the device.

The detector probe should be cleaned and disinfected after each use. Disinfection with a soft cloth soaked in ethyl or isopropyl alcohol (70–90%) is recommended.

08. TROUBLESHOOTING

If the following problems occur during use, follow the instructions below. If you cannot resolve the problem, contact an authorised service centre.

Problem	Possible cause	Solution
No sound	Battery is dead or insufficiently charged.	Replace the battery with a new one.
	Damage to the device or battery connection cable.	Check the device for visible damage. If the problem persists, contact an authorised service centre.

Quiet sound	Volume too low.	Increase the volume.
	Low battery charge.	Replace the battery.
	No or insufficient coupling gel.	Apply an appropriate amount of coupling gel to the acoustic surface of the probe.
Noise	The probe is too close to the main unit.	Keep the probe at a suitable distance from the main unit.
	Interference from external devices	Move the device away from sources of interference (e.g. mobile phones).
	Low battery level.	Replace the battery.
Low sensitivity	Incorrect probe position.	Change the position of the probe to obtain a clearer signal.
	No or insufficient conductive gel.	Apply a sufficient amount of conductive gel.

09. SPECIFICATIONS:

ENVIRONMENTAL		
Operating environment	Temperature	10°C ~+ 40°C (+50°F ~ +104°F)
	Relative humidity	15% RH ~ 85% RH
	Atmospheric pressure	70 kPa ~ 106 kPa
Transport and storage	Temperature	-20°C ~ +55°C (-4°F ~ +131°F)
	Relative humidity	10% RH ~ 95% RH
	Atmospheric pressure	58 kPa ~ 106 kPa
PHYSICAL		
Dimensions	Main unit: 136mm × 70mm × 41mm	
Weight	Main unit: 220 g (with batteries and probe)	
Screen	OLED, 2.7 × 1.5 cm	
Power	Input	3 V DC
	Operating frequency	2.0 MHz ± 5%
	Battery	2 x 1.5 V AA batteries
Type of protection against electric shock	Internal power source	
Degree of protection against electric shock	Type BF parts used	
Degree of protection against harmful water ingress	Main unit	IP22
	Probe	IPX7
Degree of safety in the presence of flammable gases	The device is not suitable for use in the presence of flammable gases	
Requirements for ultrasonic transmission medium	pH: 5.5 – 8	
OPERATING		
Operating mode	Continuous wave (CW) Doppler	
FHR – Overall sensitivity	≥90 dB	
FHR measurement range	60–210 bpm (beats per minute)	
Accuracy	±2 bpm (beats per minute)	
FHR resolution	1 bpm (beats per minute)	
Nominal ultrasound frequency	2 MHz	
Operating frequency	(2.0 ± 5%) MHz	
Effective radiation area of the probe	157 mm ² ±30%	

2 MHz wired probe	$I_{oh} < 20 \text{ mW/cm}^2$
	$p < 0.1 \text{ MPa}$
	$I_{spz} < 100 \text{ mW/cm}^2$
	$I_{sat} < 20 \text{ mW/cm}^2$
	W0<50mW
Sensitivity	From the 12th week of pregnancy
Device lifespan	3 years

10. WARRANTY CARD

Dear Customer, thank you for purchasing our Neno Sereno foetal heart rate monitor. If you encounter any problems with the device under normal conditions, please contact an authorised Neno service centre or distributor. Keep your warranty card in case repairs are necessary.

The product is covered by a 24-month warranty. The warranty terms and conditions can be found at [ronie: https://neno.pl/gwarancja](https://neno.pl/gwarancja)

Details, contact information and service address can be found at: <https://neno.pl/kontakt>

Specifications and kit contents are subject to change without notice. We apologise for any inconvenience.

APPENDIX 1: ELECTROMAGNETIC COMPATIBILITY

1. Electromagnetic emissions

Guidelines and manufacturer's declaration – Electromagnetic emissions		
The ultrasonic detector is designed to operate in the electromagnetic environment specified in the following guidelines. The purchaser or user of the device should ensure that it is used in such an electromagnetic environment.		
RF emissions CISPR 11	Group 1	The ultrasonic detector uses RF energy only for its internal functions. Therefore, RF emissions are low and should not cause interference with nearby electronic equipment.
RF emissions CISPR 11	Class B	The ultrasonic detector is intended for use in all types of premises, including residential buildings and premises directly connected to the public low-voltage power supply network supplying buildings intended for residential purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

2. Electromagnetic immunity

Guidelines and manufacturer's declaration – Electromagnetic immunity			
The ultrasonic detector is designed to operate in the electromagnetic environment specified in the following guidelines. The purchaser or user of the device should ensure that it is used in such an electromagnetic environment.			
Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidelines
Electrostatic discharge (ESD) IEC 61000-4-2	$\pm 8 \text{ kV contact}$ $\pm 2 \text{ kV}, \pm 4 \text{ kV}, \pm 8 \text{ kV}, \pm 15 \text{ kV air}$	$\pm 8 \text{ kV contact}$ $\pm 2 \text{ kV}, \pm 4 \text{ kV}, \pm 8 \text{ kV}, \pm 15 \text{ kV air}$	Floors should be wooden, concrete or ceramic tile. If floors are covered with synthetic material, relative humidity should be at least 30%.
Fast electrical surges/impulses IEC 61000-4-4	$\pm 2 \text{ kV}$ for power supply line Repetition frequency 100 kHz $\pm 1 \text{ kV}$ for input/output line	Not applied	Not applied
Surge IEC 61000-4-5	$\pm 0.5 \text{ kV}, \pm 1 \text{ kV}$ differential line-line mode	Not applicable	Not applied

Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	0% UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single-phase: at 0° 0% UT; 250/300 cycles	Not applied	Not applied
Power supply frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m, 50/60 Hz	30 A/m, 50/60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: UT is the AC supply voltage before applying the test level.

3. Electromagnetic immunity

Guidelines and manufacturer's statement – Electromagnetic immunity			
The ultrasonic detector is intended for use in the electromagnetic environment specified in the following guidelines. The purchaser or user of the device should ensure that it is used in such an electromagnetic environment.			
Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidelines
Conducted radio waves IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms 150 kHz to 80 MHz outside ISM bands	Not applicable	Portable and mobile radio communication equipment should not be used within the recommended separation distance from any part of the ultrasound detector, including its cables. This distance should be calculated using the formula appropriate for the transmitter frequency. Recommended separation distance.
Radiated radio waves IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	10V/m	$d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ <p>800 MHz to 2.7 GHz</p> <p>where P is the maximum transmitter output power in watts (W) as specified by the manufacturer, and d is the recommended separation distance in metres (m). The field strength generated by stationary RF transmitters, as determined by electromagnetic environment testing, should be less than the established compliance level in each frequency range. Interference may occur in the vicinity of devices marked with the following symbol:</p> 

NOTE: At frequencies of 80 MHz and 800 MHz, a higher frequency range applies.

NOTE: These guidelines may not apply in all situations. Electromagnetic radiation is affected by absorption and reflection by structures, objects and people.

^a Advantages of fixed transmitting stations, such as base stations for amateur (mobile/wireless) telephones and terrestrial radiotelephones, amateur radio, AM (amplitude modulation) and FM (frequency modulation) transmissions, and television transmissions, whose field strength is not theoretically possible to predict accurately. To evaluate the electromagnetic environment of a fixed RF transmitter, an electromagnetic site survey should be conducted. If the measured field strength of the foetal heart rate detector is higher than the RF compliance level of the above application, the device should be observed to see if it is functioning properly. If malfunction is observed, additional measures may be necessary, such as reorienting or repositioning the detector.

^b Over the entire frequency range from 150 kHz to 80 MHz, the field strength should be less than 3 V/m.

4. Recommended separation distances

Recommended distances between portable and mobile RF communications equipment and the device

This device is intended for use in an electromagnetic environment in which radio interference is controlled. The user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device, as recommended below, depending on the maximum output power of the communications equipment.

Maximum rated transmitter power W	Distance depending on the transmitter frequency m		
	150 kHz to 80 MHz $d = \frac{3.5}{\sqrt{P}} \sqrt{P}$	80 MHz to 800 MHz $d = \frac{3.5}{\sqrt{P}} \sqrt{P}$	800 MHz to 2.7 GHz $d = \frac{7}{\sqrt{P}} \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.20	1.2	2.3
10	3.80	3.8	7.3
100	12.00	12	23

For transmitters with a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the formula given in the relevant transmitter frequency column, where P is the maximum output power of the transmitter specified by the manufacturer in watts (W) as the unit.

Note 1: At frequencies of 80 MHz and 800 MHz, the formula for the higher frequency range is used.

Note 2: These guidelines may not be appropriate in all situations. The propagation of electromagnetic waves is affected by absorption and reflection from buildings, objects and the human body.

APPENDIX 2: INTENSITY AND SAFETY OF ULTRASOUND

5. Ultrasound in medicine

The use of ultrasound diagnostics has proven to be a valuable tool in medical practice. Given the known benefits of non-invasive medical examination and diagnosis, including human foetal examination, the question of clinical safety with regard to ultrasound intensity arises.

There is no easy answer to the question of safety associated with the use of ultrasound diagnostic equipment. The application of the ALARA (As Low As Reasonably Achievable) principle serves as a practical guideline to help achieve reasonable results at the lowest possible ultrasound power.

The American Institute of Ultrasound in Medicine (AIUM) states that, given more than 25 years of use and no confirmed biological effects on patients or device operators, the benefits of prudent use of ultrasound diagnostics clearly outweigh any risks.

6. Ultrasound safety and the ALARA principle

Ultrasound waves dissipate energy in the form of heat and can therefore cause tissue heating. Although this effect is very small with Doppler waves, it is important to know how to control and limit patient exposure. Major regulatory bodies in the field of ultrasound have issued statements that there are no known adverse effects of diagnostic ultrasound, but exposure levels should always be limited to the lowest reasonably achievable level (ALARA principle).

7. Explanation of MI/TI

a) MI (mechanical index)

Cavitation occurs when an ultrasound wave passes through and contacts tissue, causing instantaneous local heating. This phenomenon depends on the acoustic pressure, spectrum, focusing, transmission mode, and factors such as the condition and properties of the tissue and the boundary. This mechanical biological effect is a threshold phenomenon that occurs when a certain level of ultrasound power is exceeded. This threshold is related to the type of tissue. Although no confirmed adverse mechanical effects on patients or mammals have been reported due to exposure to intensities typical of current diagnostic ultrasound devices, the cavitation threshold remains undefined. In general, the higher the acoustic pressure, the greater the

potential mechanical effects; the lower the acoustic frequency, the greater the potential mechanical effects. The AIUM and NEMA have developed a mechanical index (MI) to indicate the potential for mechanical effects. The MI is defined as the ratio of peak acoustic pressure (to be calculated based on tissue acoustic attenuation of 0.3 dB/cm/MHz) to acoustic frequency:

$$MI = \frac{P_{r,a}}{f_{awt} \times C_{MI}}$$

$$C_{MI} = 1 \text{ (MPa / MHz)}$$

b) TI (thermal index)

Tissue heating is caused by the absorption of ultrasound during its application. The temperature increase depends on the acoustic intensity, the area of exposure, and the thermophysical properties of the tissue. To indicate the potential for temperature increase due to thermal effects, the AIUM and NEMA have developed a thermal index (TI). It is defined as the ratio of total acoustic power to the acoustic power required to increase tissue temperature by 1°C. Depending on the different thermophysical properties of the tissue, TI is divided into three types: TIS, TIB and TIC.

- TIS (tissue index for soft tissues): provides an estimate of the potential temperature rise in soft tissues or similar tissues.
- TIB (bone thermal index): provides an estimate of the potential temperature increase when the ultrasound beam passes through soft tissue and the focal area is in close proximity to bone.
- TIC (skull bone thermal index): provides an estimate of the potential temperature rise in the skull or superficial bones.

c) Measurement uncertainties

The measurement uncertainties were mainly systematic in nature; random uncertainties were insignificant in comparison. The overall systematic uncertainties were determined as follows:

- Hydrophone sensitivity: ±23 per cent for intensity, ±11.5 per cent for pressure. Based on the hydrophone calibration report prepared by ONDA. Uncertainty specified within ±1 dB in the 1–15 MHz frequency band.
- Digital transducer: ±3 per cent for intensity, ±1.5 per cent for pressure.
- Based on the specified accuracy of the 8-bit resolution of the Agilent DSO6012 digital oscilloscope and the signal-to-noise ratio of the measurement.
- Temperature: ±1 per cent
- Based on water bath temperature variation of ±1 °C.
- Spatial averaging: ±10 per cent for intensity, ±5 per cent for pressure.
- Non-linear distortion: not applicable.

No effects of non-linear propagation were observed.

Since all of the above sources of error are independent, they can be summed based on the effective value, resulting in a total uncertainty of ± 24.65% for all given intensity values, ± 12.33% for all pressure values, and ± 12.53% for the mechanical index.

8. Statement on prudent use

Although no confirmed biological effects on patients have been reported due to exposure to current diagnostic ultrasound devices, it is possible that such biological effects may be identified in the future. Therefore, ultrasound should be used with caution. High acoustic power levels and long exposure times should be avoided when obtaining the necessary clinical information.

9. Information on acoustic performance and safety - Bibliography

- a) "Bio-effects and Safety of Diagnostic Ultrasound" issued by AIUM in 1993
- b) "Medical Ultrasound Safety" issued by AIUM in 1994
- c) „Acoustic Output Measurement Standard for Diagnostic Ultrasound Equipment, Revision" issued by AIUM/NEMA in 2004
- d) „Standard for real-time display of thermal and mechanical acoustic output indices on diagnostic ultrasound equipment, Revision 2" issued by AIUM/NEMA in 2004
- e) „Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Probes" issued in 2008.
- f) "Medical electrical equipment—Part 2-37: Particular requirements for the basic safety and essential performance of ultrasonic medical diagnostic and monitoring equipment" issued by IEC in 2007.

10. List of acoustic probe output parameters

Table of acoustic results for track 1 Mode without automatic scanning

Operating mode: CW mode

Operating frequency: 2.0 MHz

Acoustic output		MI	$I_{\text{rms},3}$ (mW/cm ²)	$I_{\text{rms},2}$ (W/cm ²)	
Global maximum value		0.013	6.55	0.0066	
Related acoustic parameter	$P_{r,3}$ (MPa)	0.018			
	W_0 (mW)		18.22	18.22	
	f_c (MHz)	2.00	2.00	2.00	
	Z_{sp} (cm)	3.00	3.00	3.00	
	Beam dimensions	X_g (cm)		0.52	0.52
		Y_g (cm)		0.76	0.76
	PD (usec)	N/A		N/A	
	PRF (Hz)	N/A		N/A	
EBD	A_2 (cm)		Φ1.41		
	E_2 (cm)		Φ1.41		
Operational control conditions	Focus (mm)		Constant		
	Depth (mm)		Fixed		
	Frequency (MHz)		2.00		

Acoustic power reporting table for IEC60601-2-37 (IEC60601-2-37, edition 2.1, 2015-06, table 201.103)

Operating mode: CW mode

Operating frequency: 2.0 MHz

INDEX		MI		TIS		TIB		TIC
		On the surface	Below the surface	On the surface	Below the surface			
Maximum index value		0.013	0.057		0.13			N/A
Index component value			N/A	0.057	N/A	0.13		
Acoustic parameters	$p_{r,a}$ at $z_{p,r,a}$ (MPa)	0.018						
	P (mW)		18.22		18.22			N/A
	$P_{r,3}$ (mW)			N/A	N/A			
	z (cm)			3.0				
	z_b (cm)					3.0		
	$z_{p,r}$ (cm)	3.0						
	$z_{p,r,a}$ (cm)	3.0						
	$f_{z,ref}$ (MHz)	2.0	2.0		2.0			N/A
Other information	p_{rr} (Hz)	N/A						
	s_{rr} (Hz)	N/A						
	$n_{p,r}$	N/A						
	$I_{p,r,a}$ at $z_{p,r,a}$ (W/cm ²)	0.0066						
	$I_{sp,r,a}$ at $z_{p,r,a}$ or $z_{p,r}$ (mW/cm ²)	6.55						
	$I_{sp,r,a}$ at $z_{p,r}$ or $z_{p,r}$ (mW/cm ²)	9.90						
	p_r at $z_{p,r}$ (MPa)	0.022						
	Operating control conditions	Sharpness (mm)	Constant					
Depth (mm)		Constant						
Frequency (MHz)		2.00						

APPENDIX 3 OVERALL SENSITIVITY

Overall sensitivity (2 MHz wired probe)

Target reflector diameter (mm)	Distance (d) (mm)	Reflection loss A(d)	Bidirectional attenuation $B = \Sigma B_s + B_w$			V_s (r.m.s) mV	V_n (r.m.s) mV	Overall sensitivity (S=A(d)+B+C) dB				
			ΣB (T: ultrasonic attenuation model No., dB)	B_s (dB)	B (dB)							
1.58 A=45.7dB@ 2MHz	50	39.8	T	6#6#6#1#	0	78.3	70.06	31.43	6.96	125.06		
			B _s	78.3								
	75	43.3	T	6#6#6#	0	74.6	69.39	29.83	7.33		125.23	
			B _s	74.6								
	100	45.8	T	6#6#5#4#	0	72.6	69.72	29.66	7.42			125.82
			B _s	72.6								
200	51.8	T	6#6#4#3#	0	67.6	54.04	29.35	5.30	124.70			
		B _s	67.6									
50	37.2	T	6#6#6#2#	0	80.0	60.47	31.09	5.78		122.98		
		B _s	80.0									
75	40.7	T	6#6#6#1#	0	78.3	58.64	30.32	5.73			124.73	
		B _s	78.3									
100	43.2	T	6#6#6#1#	0	78.3	63.61	30.74	6.32	127.82			
		B _s	78.3									
200	49.2	T	6#6#4#3#	0	67.6	52.58	29.31	5.08		121.88		
		B _s	67.6									
Detector frequency (Hz)		200			Target speed (cm/s)		10					

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MedNet EC-REP GmbH
Shenzhen IMDK Medical Technology CO.Ltd