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




Medic T02

A**B****C****D****E****F****G**

H

	1
	2
	3
	4
	5
	6
	7
	8
	9
	10
	11
	12
LOT	13
	14
	15
IP22	16
	17
CE 0123	18
	19
MD	20
	21

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Szanowny Kliencie,

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Dear customer,
thank you for purchasing the NENO MEDIC T02 AET-R161. The device you have received is a non-contact, reusable infrared thermometer. The device can be used both at home and by professionals. The thermometer measures body temperature in adults and children over 3 months of age. Before using the device, read the instructions below. Keep these instructions for later reference.

01. INTRODUCTION

1. Try the device on yourself to understand how it works.
2. To avoid affecting the accuracy of the device due to an external environment that is too cold or too hot, it should be placed in a general environment with a room temperature of 15°C-40°C for at least 30 minutes before use to avoid incorrect measurement results.
3. The person whose temperature is being measured should be in a room at room temperature for at least 20 minutes before the measurement is taken.
4. After exercise, wait 30 minutes before taking a temperature measurement.
5. Wipe the skin whose temperature is being measured and ensure that the hair does not obstruct the sensor on the subject's forehead.
6. If you take multiple measurements in a row, wait 5 seconds between measurements and take them from different distances to achieve the most accurate measurements.
7. Ensure that the sensor is not contaminated before each measurement.
8. The standard body temperature, in both adults and children, is around 36.6 degrees Celsius. Consult your doctor if the measured temperature is significantly higher.

Intended use/Indications for use

The non-contact thermometer is designed to measure the body temperature of people aged 3 months and above, detecting infrared heat during non-contact measurement (approximately 1.5-5.0 cm) from the forehead in a home or hospital environment.

02. METHOD OF APPLICATION

SEE FIG. A

The thermometer measures the temperature of a person's forehead and works from 1.5 cm to 5 cm.

1. Switch on the thermometer by pressing the **SEE FIG. E** button. When you hear a beep and the screen displays ---°C, the unit is ready for use.
2. Place the thermometer at a suitable distance from the forehead of the person to be measured, with the sensor pointing towards the forehead. Press and wait a second until you hear a short beep. The display will show the temperature measurement taken. If you want to perform multiple tests in a row, wait 5 seconds between measurements. This procedure can be repeated several times.
3. To switch off the unit, hold **SEE FIG. E** for approximately 5 seconds. The thermometer will automatically switch off after 60s+/-10s.
4. If the measurement is below 32°C, the display will show „Lo” and the thermometer will beep twice.
5. If the measurement is above 37.8°C, the thermometer will give off six signals.
6. If the measurement is above 42.2°C, the display will show „Hi” and the thermometer will beep twice.

03. CHANGE OF UNITS OF MEASUREMENT

To change the unit of measurement from degrees Celsius to degrees Fahrenheit, hold the **SEE FIG. E** button for 8 seconds while the thermometer is off until the screen displays ---°C. Press **SEE FIG. E** to change the unit of measurement. The thermometer will automatically switch off after 4 seconds.

04. MEASUREMENTS STORED IN MEMORY

1. Make sure the thermometer is switched off.
2. Hold down **SEE FIG. E** for 4 seconds to activate memory mode. The last measurement taken by the thermometer will be displayed.
3. Press the **SEE FIG. E** button to toggle between the measurements stored in memory. The higher the number displayed before the temperature, the older the measurement displayed.
4. The thermometer can show the last 32 measurements taken.
5. In memory mode, the thermometer switches off spontaneously after 12 seconds of inactivity.

05. DESCRIPTION OF THE DEVICE

SEE FIG. B

1. Sensor
2. Housing
3. Space for batteries
4. Display
5. On/off button

06. DESCRIPTION OF THE DISPLAY FUNCTIONS

SEE FIG. C

1. MEASUREMENT PLACE - indicates the location from which the temperature measurement should be taken (patient's forehead).
2. UNIT OF TEMPERATURE - shows whether the temperature is measured in degrees Celsius or Fahrenheit.
3. MEASURED TEMPERATURE - shows the temperature measurement after it has been taken.
4. STORED IN DEVICE MEMORY INDICATOR - shows that the measurement currently being viewed is an archive measurement, stored in the device's memory.
5. BATTERY INDICATOR - Appears when the batteries are close to running low and reminds you to replace them.

07. BATTERY REPLACEMENT

SEE FIG. D

1. When the battery indicator icon appears on the display, immediately replace the batteries with two new ones.
2. Press the lid of the battery compartment with your thumb and slide it open to create a gap between the compartment and the housing of the device. Pull the lid up and to the side to remove it.
3. Remove the old batteries and insert two new AAA batteries. Pay attention to the poles of the batteries so that they are not placed in reverse.
4. Place the battery compartment cover back in place.

Note: Remove the batteries from the device if you will not be using it for a long period of time. Dispose of batteries only in the containers intended for them.

08. PRECAUTIONS

SEE FIG. G

1. Do not use the thermometer in rooms where the temperature is below 15°C(59°F) or above 40°C(104°F) or where the humidity is below 10%RH or above 80%RH. These conditions may cause measurement errors.
2. Do not expose the thermometer to temperatures below -20°C or above 55°C or humidity below 15%RH or above 95%RH. These conditions may cause damage to the device.
3. If the device's sensor is damaged or loose, contact service immediately.
4. Do not touch the sensor or blow on it. This can lead to incorrect

measurements.

5. Do not self-diagnose using the measurements obtained. If the results obtained concern you, consult your doctor.
6. If the sensor is dirty, gently wipe it with a soft material. Cleaning the sensor with paper towels or toilet paper can scratch and damage it.
7. Incorrect installation of the battery in the device will lead to excessive heating of the battery and damage to the device
8. Do not expose the device to liquids. The device is not waterproof.
9. The sensor is made of a special type of glass which is difficult to damage. If you notice any damage, stop using the device immediately and contact your dealer.
10. Do not use a mobile phone near a switched-on thermometer. Portable communication equipment can affect the measurement of medical electronic equipment.
11. Do not modify the appliance yourself or take it apart. This will void the warranty.
12. Keep the thermometer out of the reach of children. Some parts of the device are so small that a child could swallow them. If a child accidentally swallows the battery or protective film of the device, consult a doctor immediately.
13. If the appliance has been modified, carry out the appropriate tests to ensure that it is functioning correctly. If problems occur, consult service immediately.
14. If there is a temperature difference between the room where the thermometer is stored and the room where the measurement is taken, make sure that the thermometer is kept in the measurement room for at least 30 minutes before the first measurement is taken. This will avoid anomalies.
15. The device does not require periodic calibration. In the event of measurement errors, contact your dealer.
16. The person using the thermometer should not touch the thermometer and the patient at the same time. This can lead to measurement errors.
17. Do not repair or clean the appliance yourself while it is in use.

The patient can self-measure and replace the batteries by following the guidelines in this manual.

09. MAINTENANCE, STORAGE AND CALIBRATION OF THE INSTRUMENT

Maintenance:

1. Remove stains from the unit with a soft, dry material.
2. Gently wipe off dirt on the device with a fine material, use a cotton swab or a fine material with alcohol to gently wipe off the sensor and display.
After use, the product can be sterilised using a piece of cotton or a soft cloth soaked in medical alcohol (70-80% alcohol content).
3. Do not wash the unit with water or cleaners containing caustic agents or benzene.
4. Do not immerse in any liquids.

Storage:

1. The device must not be stored where: the temperature is too high or too low, the humidity is too high or too low, the device is exposed directly to sunlight, near the flow of electricity or in dusty areas. These conditions may damage the device and cause inaccurate measurements.
2. Remove the batteries from the device if it will not be used for a long time.

Calibration:

1. The thermometer is calibrated at the time of manufacture. If the thermometer is used in accordance with the instructions, it does not require periodic calibration. If you doubt the accuracy of the measurement, please contact your dealer or service centre.
Do not attempt to modify or disassemble and fold the thermometer yourself

10. SPECIFICATION

Device name: Infrared thermometer | **Measurement mode:** forehead | **Units of measurement:** degrees Celsius (°C) or Fahrenheit (°F) | **Opera-**

ting conditions: Temperature: 15°C-40°C; **Relative humidity:** 10%-80%; **Atmospheric pressure:** 70-106Kpa | **Storage and transport conditions: Temperature:**-20°C-55°C; **Relative humidity:** 15%-93%RH; | **Atmospheric pressure:** 70-106Kpa | **Thermometer time from highest or lowest environmental temperature to operating environment:** 30 minutes | **High temperature indication:** ≥37.8°C (100°F) | **Software version:** A.01 | **Measurement distance:** 15 to 50mm from the measurement point | **Measurement range (Frontal mode):** 32.0°C to 42.2°C(89.6°F-108°F) | **Measurement accuracy (for laboratory conditions):** ±0.2°C (35.0°C-42.0°C)/±0.4°F (95°F-107.6°F), ±0.3°C (±0.5°F) remaining range | **Clinical repeatability:** ≤±0.3°C (0.5°F) | **Accuracy for clinical test:** ≤±0.3°C (0.5°F) | **Display resolution:** 0.1°C (0.1°F) | **Measurement time:** 1s (5s later for the next measurement) | **Battery life:** New batteries (2 AAA batteries) will take approx. 1000 measurements | **Auto shut-off function:** 60s±10s | **Memory:** 32 last measurements taken | **Dimensions:** 150x37x16.3mm | **Weight:** 51g (without batteries) | **Batteries:** 2xAAA (DC 3V) | **Product lifetime:** 5 years | **Measurement location:** forehead | **Reference side:** Mouth

Note: Not intended for sterilisation.

Do not use in oxygen-rich environments.

It is not intended for use in an emergency medical services environment.

It is not intended for use with flammable agents.

It is not intended for use with flammable anaesthetics.

11. KIT CONTENTS:

Non-contact thermometer x1

User manual x1

12. EXPLANATION OF SYMBOLS

The following symbols may appear in the instructions, on the thermometer or on the rest of the kit.

SEE FIG. H

1. Indicates information with which it is mandatory to be familiar before using the device.
2. Follow the instructions for use.
3. BF type part application.
4. Note: Read the accompanying documentation.
5. Do not dispose of the product in the container for mixed municipal waste. Dispose of the product in accordance with the guidelines for the disposal of electronic devices of this type.
6. The transport packaging should be kept out of the rain.
7. The transport packaging should not be exposed to sunlight.
8. Indicates the side of the transport packaging that should be facing upwards.
9. The contents of the transport packaging are fragile, the packaging should be handled with special care.
10. Indicates the temperature range at which the transport packaging should be stored.
11. Non-ionising electromagnetic radiation.
12. Don't roll over
13. LOT number
14. Date of production
15. The product cannot be used after the date listed next to this icon.
16. The product is resistant to solid objects with a diameter of 12.5mm or more and against vertically falling water droplets when the unit is raised at an angle of up to 15 degrees.
17. Manufacturer
18. CE mark: The product complies with EU 2017/745
19. Authorised representative in the European Community.
20. Medical device.
21. Template UDI
 - (01) DI Identification of the device
 - (17) Expiry date
 - (11) Production date
 - (10) Lot number

13. ELECTROMAGNETIC COMPATIBILITY

MEDICAL EQUIPMENT is suitable for use for healthcare purposes in home and hospital environments.

Caution: Do not operate the device in the vicinity of high-frequency surgical equipment or in rooms designed for electromagnetic resonance (with electromagnetic shielding), as these areas are subject to high levels of interference electromagnetic.

Caution: Using the device in close proximity to other devices can lead to incorrect measurements and should be avoided. If the use of two devices in close proximity to each other is unavoidable, both devices should be closely monitored to ensure they are operating correctly.

If present: list of cables and maximum cable lengths (if applicable), transducers and other ACCESSORIES that may be specified by the RESPONSIBLE ORGANISATION and that may affect the compliance of the MEDICAL EQUIPMENT or the MEDICAL SYSTEM with the requirements of Clause 7 (EMISSIONS) and Clause 8 (RESISTANCE). ACCESSORIES may be specified in general (e.g. shielded cable) or specifically (e.g. specifying the PRODUCER and TYPE OF EQUIPMENT).

If it is: The functioning of the MEDICAL EQUIPMENT or MEDICAL SYSTEM has been identified as KEY FUNCTIONING and the changes that the USER may experience if the KEY FUNCTIONING is disrupted by ELECTROMAGNETIC INJURY are described (the term KEY FUNCTIONING need not be used).

1. All necessary instructions for maintaining the BASIC SAFETY and KEY FUNCTION of the equipment in relation to electromagnetic interference during the serviceability period.
2. Manufacturer's information and declaration - electromagnetic emissions and immunity to disturbance.

Table 1

Guidance and manufacturer's declaration – electromagnetic emission	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	N/A
Voltage fluctuations/flicker emissions IEC 61000-3-3	N/A

Table 2

Guidance and manufacturer's declaration - electromagnetic Immunity		
Immunity test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	Power supply lines : ±2 kV input/output lines : ±1 kV	N/A
Surge IEC 61000-4-5	line(s) to line(s) : ±1 kV. line(s) to earth : ±2 kV. 100 kHz repetition frequency	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% 1 cycle And 70% 25/30 cycles Single phase: at 0 0% 300 cycle	N/A
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz

Conducted RF IEC61000-4-6	150KHz to 80MHz : 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	N/A
Radiated RF IEC61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m 80MHz - 2.7GHz 80% AM at 1kHz
NOTE: Ut is the AC mains voltage before applying the test level.		

Table 3

Information and manufacturer's declaration - resistance to disturbance							
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	Distance (m)	Immunity Test level (v/m)
	385	380-390	TETRA 400	Pulse modulation 18Hz	1,8	0,3	27
	450	380-390	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0,3	28
	710 745 780	704-787	LTE Band 13, 17	Pulse modulation 217Hz	0,2	0,3	9
	810 870 930	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18Hz	2	0,3	28
	1720 1845 1970	1700-1990	GSM 1800; C D M A 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0,3	28
	2450	2400-2570	Bluetooth, WLAN, 802.11 b/ g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0,3	28
	5240 5500 5785	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0,3	9

14. ERROR MESSAGES

Communication	Problem	Solution
Hi	Temperature above 42.2°C	Use the thermometer only under conditions specified in the manual. If the error persists, consult the service.
Lo	Temperature below 32°C	Use the thermometer only under conditions specified in the manual. If an error persists, consult the service.
Err	The sensor is not fixed	Contact service
ErE	Memory EEPROM no operates correctly	Contact service
SEE FIG. F	Discharged batteries	Replace new batteries
ErH	Ambient temperature too high	Decrease the temperature of the environment, keep it within 15-40°C
ErL	Ambient temperature too low	Decrease the temperature of the environment, keep it within 15-40°C

15. WARRANTY CARD

The product is covered by a 24-month warranty. The warranty conditions can be found at: <https://neno.pl/gwarancja> Details, contact and service address can be found at: <https://neno.pl/kontakt>

Specifications and contents are subject to change without notice. We apologise for any inconvenience. KGK Trend declares that the NENO Medic T02 complies with the essential requirements of Directive 2014/53/EU. The full declaration can be found at the following link: <https://neno.pl/download/DOC/deklaracja-CE-Neno-Medic-T02.pdf>

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EN

The crossed out trash can symbol indicates that unusable electrical or electronic devices, its accessories (such as power supplies, cords) or components (for example batteries, if included) cannot be disposed of alongside with household waste. In order to dispose of the devices or its components (for example, batteries) deliver the device to the collection point, where it will be accepted free of charge. Disposal is subject to the recast version of the WEEE Directive (2012/19/ EU) and the Directive on batteries and accumulators (2006/66 / EC). Proper disposal of the device prevents degradation of the natural environment. Information about the collection points of the facilities is issued by the competent local authorities. Incorrect disposal of waste is subject to penalties provided for by the law in force in the given area.

neno[®]

Importer:

KGK TREND Sp. z o.o.
st. Ujastek 5b, 31-752 Cracow, Poland
Made in P.R.C.

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Company name: CMC Medical Devices & Drugs S.L.
Address: C/ Horacio Lengo Nº 18, CP 29006, Málaga, Spain
Web: www.cmcmedicaldevices.com
Tel: +34951214054 Email: Info@cmcmedicaldevices.com

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ALICN MEDICAL SHENZHEN, INC
Room 410, Building A, 3rd Sub-park, Leibo Zhongcheng Life Science Park,
No. 22 Jinxiu East Road, Jinsha Community, Kengzi Street,
Pingshan District, 518118 Shenzhen, Guangdong,
PEOPLE'S REPUBLIC OF CHINA

The quality of the product is guaranteed by the TÜV SÜD Product Service GmbH

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Jakość produktu gwarantuje TÜV SÜD Product Service GmbH