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


Medic T07

| | |
|--|-----|
|  | A.1 |
|  | A.2 |
| IP22 | A.3 |
|  0598 | A.4 |
|  | A.5 |
|  | A.6 |
|  | A.7 |
| EC REP | A.8 |
| MD | A.9 |

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USER MANUAL

Dear Customer,

Thank you for choosing Neno Medic T07. The device you have purchased is a non-contact thermometer that measures body or object temperature using an infrared light wave sensor. Please read the instructions below before use.

01. PRECAUTIONS

- Do not use the non-contact thermometer for any purpose other than that described in the instructions. The thermometer is suitable for both domestic use and as a medical product.
- Do not immerse the thermometer in water or other liquids. When cleaning the device, refer to the instructions under „Cleaning and Storage“.
- The thermometer should be stored in a dry, clean place, away from the sun. The thermometer works best at a temperature of 15-40°C and at a humidity of 15-95%RH.
- Do not touch the thermometer sensor.
- Sweat, hair, headgear etc. may underestimate the measured temperature. Make sure that nothing obstructs the sensor on the subject's bare skin.
- Do not drop the product, take it apart or carry out repairs or modifications yourself.
- Do not keep the thermometer near strong electrostatic fields or magnetic fields that may cause measurement errors.
- If problems occur, stop using the device and contact your dealer.
- Do not dispose of this product or its batteries in your household waste. Follow the laws that apply to the disposal of electronic equipment and batteries.
- If the device will not be used for a long time, remove the batteries to avoid the risk of damaging the thermometer.
- Do not put new and partially used batteries in the device at the same time. This can damage the device.
- The infrared thermometers take human body temperature via the eardrum or forehead for both professional and home use. Intended patient population for forehead mode applies to all age group and ear mode applies to people above 3 months.

Caution! Keep thermometer out of reach of children. Do not dispose of batteries in a fire. The thermometer is not a substitute for medical examination and recommendations.

02. EXPLANATION OF SYMBOLS

SEE FIGURES A.1-A.8

A.1 Device with BF-type parts | **A.2** Do not dispose of the product in the mixed municipal waste container. Dispose of the product in accordance with the guidelines for the disposal of electronic devices of this type | **A.3** Luminaires with condensation and drip-proof protection | **A.4** CE mark Product meets EU requirements | **A.5** Follow instructions for use | **A.6** Manufacturer | **A.7** Date of manufacture | **A.8** Authorised representative in the European Community | **A.9** Medical device

03. PRODUCT DESCRIPTION

Purpose of the device

The non-contact thermometer is used to measure body temperature (forehead, ear) and object temperature. The device is suitable for both home use and for use as a medical device. The thermometer can be used to measure temperature regardless of the age of the subject. Do not use the thermometer if the ear is infected with otitis or suppuration. The device has no side-effects if administered correctly and residual risk is acceptable.

Construction of the device: **SEE FIG. B**

1. LED display
2. Power and measurement button
3. Memory/Mute button - press the button to view stored measurements/hold the button until the crossed-out speaker icon appears or disappears to mute the device or restore sounds
4. Mode button - pressing this button toggles the device between body and object temperature mode and between adult and child modes up to 12 years of age
5. Probe
6. Button to change the unit (°C/°F)
7. Probe cover - when the probe cover is removed the thermometer automatically switches to ear measurement mode. Similarly, when the probe cover is placed on the device, the body/object temperature will be measured
8. Battery cover

Display: **SEE FIG. C**

1. Object temperature mode
2. Forehead temperature mode
3. Child mode
4. Ear temperature mode
5. Mute
6. Recall of measurements from memory
7. Temperature unit (°C)
8. Temperature unit (°F)
9. Low battery
10. Temperature value

04. USE OF THE DEVICE

1. Installing batteries:

- a. Remove the battery cover.
- b. Insert two AAA batteries. Make sure the battery terminals are facing the right way.
- c. Replace the battery cover and close the case.
NOTE: Incorrectly installed batteries can damage the thermometer.
If the batteries or device have signs of leakage or mould on them, stop using them immediately.
Do not keep batteries near a fire or throw them into a fire. This could cause an explosion.
Do not store batteries in rooms with high temperatures and moisture.
To avoid short circuits, do not keep batteries or metal objects (such as coins or keys) near electrical appliances.

2. Preparation for measurement:

- Follow the instructions below to ensure the most accurate temperature measurement:
- a. Before taking a measurement, push the hair off your forehead and clean your skin of sweat.
 - b. Select child or adult mode using the „Mode“ button.
 - c. When taking the measurement, aim the thermometer at the center of the subject's forehead, above the eyebrows. Hold the thermometer at a distance of 1~3 cm from the subject. When you press and hold the measurement button for 1 second, the temperature measurement value will be shown on the display.
 - d. If the body temperature of the person differs significantly from the temperature in the measurement room, the person should wait at least 5 minutes in the measurement room before taking the measurement.
 - e. A cold compress or other forehead cooling methods in people with a fever will mean that the measured temperature may be lower.
 - f. The temperature in the room where the measurement is taken should be stable. Do not take the measurement in rooms with high airflows such as rooms cooled by fans or ventilation systems.

- g. The thermometer should be in the same room as the measurement. If the thermometer was brought in from another room, leave it in the measurement room for at least 20 minutes before taking the measurement.
- h. Do not expose the thermometer to strong sunlight.

3. Making a measurement:

a. Measurement of body temperature:

Make sure the probe cover is on. Select forehead mode using the „Mode“ button. Aim the thermometer sensor at the forehead of the person whose temperature you want to measure. Press the measurement button for 1 second.

b. Ear temperature measurement:

Remove the probe cover from the thermometer before measuring. The thermometer will automatically change to the measurement mode when the probe cover is removed. You can use the „Mode“ button to change to baby mode. Insert the probe into the ear canal. Correct placement of the probe is essential to obtain an accurate measurement. Children under 1 year Pull the ear straight back. Children 1 year and older to adult: Pull the ear up and back. **See FIGURE D** for correct placement and take the measurement.

Note: Do not force the thermometer into the ear canal. Otherwise the ear canal may be damaged. When taking the temperature of an adult, gently pull the ear up and back to make sure that the ear canal is straight so that the temperature probe can receive infrared radiation from the eardrum. Be careful when taking the temperature of a child whose ear canal is small. Make sure that the ear canal is clean and dry before starting the measurement. In case of dirt, it is recommended to clean the ear canal. Otherwise the temperature probe may be contaminated and temperature readings may be inaccurate.

c. Measurement of the subject's room temperature:

Make sure the probe cover is on. Select the object mode using the „Mode“ button. Aim thermometer sensor at object. Press the measurement button for 1 second.

4. Viewing of stored measurements

When the unit is switched on, press the „Memory/ Mute“ button to view the measurements stored in memory. Measurement 01 is always the last measurement taken by the device. If there is no stored measurement in the device's memory, the sequence number will be displayed normally, but instead of temperature measurement, the display will show „---“. The thermometer can store up to 20 temperature measurements. If you take more measurements the oldest stored result will be deleted. Measurements taken for objects are not stored in memory.

05. ERROR MESSAGES

1. **HI** - measured temperature too high, outside the measuring range
2. **Lo** - too low a measured temperature, outside the measuring range
3. **Er1** - operating temperature outside the range 10~40°C
4. **ErC** - error occurs when data is read or written into the memory or the temperature correction is not completed
5. When the battery voltage drops below $2.5V \pm 0.1V$, the low battery symbol appears in the display. Replace the batteries.
6. Blank screen:
 - a. The thermometer switches off automatically after a short period of inactivity, this is a deliberate action of the device - restart the device by pressing the power/measurement button
 - b. Battery incorrectly installed - ensure battery terminals are facing the correct way
 - c. Batteries flat - replace with new batteries
 - d. Display does not start despite above instructions - contact service

06. CLEANING AND STORAGE

1. Do not touch or press the device sensor.
2. To clean the thermometer remove the batteries and then clean the ear-tip with a soft cloth avoiding the temperature sensor lens. The lens itself is cleaned with a cotton swab and the casing with a slightly damp, soft cloth.
3. **Note:** Keep water away from the lens during the cleaning process. Otherwise, the lens may be damaged. The lens may be scratched if cleaned with a hard object, which may cause inaccurate readings. Do not clean the thermometer with corrosive cleaning agents. Do not immerse any part of the thermometer in liquid or allow liquid to penetrate the thermometer during the cleaning process.

07. TYPICAL HUMAN BODY TEMPERATURE

The human body is a complex biological system and the range of temperature that can be considered „normal“ depends in large part on what part of the body we measure and on factors such as age, sex, skin colour and skin thickness. Women's body temperature is higher than men's by about 0.3°C. In addition, body temperature in women increases by an additional 0.3-0.5°C during ovulation

08. SPECIFICATION

Place of measurement: forehead, ear, room, food, object | **Units of measurement:** degrees Celsius (°C) or degrees Fahrenheit (°F) | **Operating temperature:** 10~40°C | **Storage temperature:** -20~50°C | **Measurement distance:** 1~3 cm | **Measurement**

range: forehead 22~43°C (71.6~109.4°F), ear 34~43°C (93.2~109.4°F), object 0~100°C (32~212°F) | **Measurement accuracy:** Forehead: ±0.2°C in the range 36.0°C-39.0°C, ±0.3°C in the ranges: 22.0°C-36.0°C and 39.0°C-43.0°C; Ear: ±0.2°C in ranges 36.0°C-39.0°C, ±0.3°C in ranges: 34.0°C-36.0°C and 39.0°C-43.0°C; Subject/ear: ±1.0°C/±2.0°F | **Memory:** 20 measurements | **Dimensions:** 145.5x45.4x43.4 mm | **Weight (without batteries):** 72g | **Batteries:** 2xAaa, DC 3V (included) | **Auto shut-off:** after 10 seconds of inactivity

09. WARRANTY CARD

Dear customer, thank you for purchasing our Neno Medic T07 thermometer. If you have any problems operating the device under normal conditions, please contact an authorized Neno brand service center or distributor. Keep your warranty card in case of repair.

Product comes with a 24-month warranty. Warranty conditions can be found at: <https://neno.pl/gwarancja>

Details, contact and service address can be found at: <https://neno.pl/kontakt>

Specifications and kit contents are subject to change without notice. We apologise for any inconvenience.

KGK Trend declares that the Neno Medic T07 device complies with the essential requirements of Directive 2014/53/EU. The full text of the declaration can be found at the link: <https://neno.pl/download/DOC/deklaracja-CE-Neno-Medic-T07.pdf>

10. ELECTROMAGNETIC COMPATIBILITY

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the ME equipment, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Table 1

| declaration - electromagnetic emission | |
|---|----------------|
| Emissions test | Compliance |
| RF emissions CISPR 11 | Group 1 |
| RF emissions CISPR 11 | Class B |
| Harmonic emissions IEC 61000-3-2 | Not applicable |
| Voltage fluctuations/ flicker emissions IEC 61000-3-3 | Not applicable |

Table 2

| declaration - electromagnetic immunity | | |
|--|--|--|
| Immunity test | IEC 60601 test level | Compliance level |
| Electrostatic discharge (ESD) IEC 61000-4-2 | ±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air | ±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air |
| Electrical fast transient/burst IEC 61000-4-4 | ± 2 kV for power supply lines ± 1 kV for input/output lines | Not applicable |
| Surge IEC 61000-4-5 | ± 0.5kV, ± 1 kV line(s) to lines ± 0.5kV, ± 1 kV, ± 2 kV line(s) to earth | Not applicable |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | 0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycles | Not applicable |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 30 A/m | 30 A/m |

NOTE: UT is the a.c. mains voltage prior to application of the test level.

Table 3

| declaration - electromagnetic immunity | | |
|--|---|------------------|
| Immunity test | IEC 60601 test level | Compliance level |
| Conducted RF IEC 61000-4-6 | 3 V 0.15 MHz to 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz | Not applicable |
| Radiated RF IEC 61000-4-3 | 10V/m 80 MHz to 2.7 GHz | 10V/m |

Table 4

| declaration - IMMUNITY to proximity fields from RF wireless communications equipment | | | | | |
|--|----------------------------------|-------------------------------|---------------|----------------|------------------|
| Immunity test | IEC60601 test level | | | | Compliance level |
| | Test frequency | Modulation | Maximum power | Immunity level | |
| Radiated RF IEC 61000-4-3 | 385 MHz | **Pulse Modulation: 18Hz | 1.8W | 27 V/m | 27 V/m |
| | 450 MHz | *FM+ 5Hz deviation: 1kHz sine | 2 W | 28 V/m | 28 V/m |
| | 710 MHz 745 MHz 780 MHz | **Pulse Modulation: 217Hz | 0.2 W | 9 V/m | 9 V/m |
| | 810 MHz 870 MHz 930 MHz | **Pulse Modulation: 18Hz | 2 W | 28 V/m | 28 V/m |
| | 1720 MHz 1845 MHz 1970 MHz | **Pulse Modulation: 217Hz | 2 W | 28 V/m | 28 V/m |
| | 2450 MHz | **Pulse Modulation: 217Hz | 2 W | 28 V/m | 28 V/m |
| | 5240 MHz 5500 MHz 5785 MHz | **Pulse Modulation: 217Hz | 0.2 W | 9 V/m | 9 V/m |

Note* - As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

Note** - The carrier shall be modulated using a 50 % duty cycle square wave signal.

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EN

The crossed out trash can symbol indicates that unusable electrical or electronic devices, its accessories (such as power supplies, cords) or components (for example batteries, if included) cannot be disposed of alongside with household waste. In order to dispose of the devices or its components (for example, batteries) deliver the device to the collection point, where it will be accepted free of charge. Disposal is subject to the recast version of the WEEE Directive (2012/19/ EU) and the Directive on batteries and accumulators (2006/66 / EC). Proper disposal of the device prevents degradation of the natural environment. Information about the collection points of the facilities is issued by the competent local authorities. Incorrect disposal of waste is subject to penalties provided for by the law in force in the given area.

Model: JPD-FR301

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Shenzhen Jumper Medical Equipment Co., Ltd.
Address: D Building, No. 71, Xintian Road, Fuyong Street,
Baoan, Shenzhen, Guangdong, China
Tel:+86-755-26696279
Fax:+86-755-26852025
Website: <http://www.jumper-medical.com>

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Importer:

KGK TREND Sp. z o. o.
Ujastek 5b, 31-752 Cracow, Poland.
Made in PRC

CE 0598

Version: 1.2
Date: 05.2024

neno[®]

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MedPath GmbH