

neno[®] beauty

ALGO

ADVANCED FACIAL FIRMING
AND CONTOURING DEVICE
WITH EMS AND RF TECHNOLOGY

Exclusive Series



EMS & RF technology

Exclusive Series

Discover the Neno Beauty series

created to support your daily
skincare. Without leaving home,
you can enjoy professional
treatments that restore your skin's
glow and freshness day by day.
With Neno Beauty, you reveal the
beauty that has always been
within you.



Exclusive Series



2 massage heads



Memory function



EMS mode




2 wavelengths of light



EMS+RF mode



3 programs



RF mode



Timer



Exclusive Series

Meet Neno Beauty Algo

Neno Beauty Algo is an advanced home skincare device that combines EMS and RF technologies, sonic vibrations, LED light, and precisely designed massage heads. Its purpose is to naturally firm and sculpt facial contours while supporting the skin's regenerative processes.



Exclusive Series




Your skincare, your ritual

You can choose from three modes – EMS, RF with yellow LED light, and EMS+RF with red LED light – allowing you to easily tailor care to your skin's needs and desired results. Each mode also offers two intensity levels. The device comes with two attachments: diamond massage heads and a flat head with eight massage points, ensuring care is fully adapted to your skin.

neno
beauty

Exclusive Series



EMS technology

Neno Beauty Algo uses EMS (Electrical Muscle Stimulation), which sends gentle electrical impulses that trigger short, controlled muscle contractions. This acts like a subtle workout for the face – helping maintain natural muscle tone and supporting a firm, youthful look. You may feel a light tingling during treatment. Two intensity levels let you adjust the device to your needs and comfort. Regular short sessions bring the best results, depending on consistency and individual skin predispositions.




Exclusive Series

A woman with dark hair and a white top is holding a white and gold RF device against her cheek. The device has a gold-colored tip and a silver band.

RF technology


Neno Beauty Algo uses radiofrequency at 4 MHz. RF waves penetrate the skin, creating gentle heat that stimulates natural regeneration and boosts collagen and elastin production – proteins essential for skin firmness and elasticity. The treatment leaves the surface intact and feels like pleasant, even warmth.



Exclusive Series

LED light therapy

Neno Beauty Algo uses carefully selected LED wavelengths. Red LED (630 nm) supports natural regeneration, making skin firmer and more elastic, while yellow LED (590 nm) has a soothing effect, helps even out skin tone, and gives a fresh, radiant look.



Exclusive Series

Two wavelengths of light

590 nm

YELLOW LIGHT

Yellow light with a wavelength of 590 nm has a soothing effect, helps even out skin tone, and gives the skin a fresh, radiant appearance.

630 nm

RED LIGHT

Red light with a wavelength of 630 nm stimulates microcirculation, helps even out skin tone and reduce discoloration, and also has a soothing and calming effect on the skin.

Exclusive Series



Two massage heads

The device includes two different attachments so you can tailor care to your skin's needs. One side features diamond massage heads, the other a flat head with eight massage points.

Diamond massage heads glide smoothly along the jawline and cheeks. Combined with EMS, they gently stimulate facial muscles, helping define the oval and highlight the V-contour exactly where you need it.

Flat head with eight massage points emits red and yellow LED light and uses sonic vibrations to stimulate microcirculation and improve skin firmness.

Alternate the attachments for a varied routine – keeping care flexible and adapted to your preferences.




Exclusive Series

EMS + diamond head mode

The first mode combines EMS with diamond head massage, offering an effective way to firm and sculpt the face. Gentle impulses plus mechanical massage strengthen muscle tone, highlight natural contours, and add skin elasticity.



Exclusive Series



EMS + RF + Red LED mode

This intensive regeneration and firming mode combines EMS with RF waves, red LED light (633 nm), sonic vibrations, and the flat head with eight massage points. It strengthens facial muscles, stimulates collagen production, and improves skin tone.


Exclusive Series



A close-up photograph of a woman with dark hair and closed eyes, using a handheld white and silver device on her cheek. The device has a glowing yellow light at its tip. The background is softly blurred.

RF + yellow LED mode


The final mode uses RF waves with yellow LED light (590 nm), sonic vibrations, and the flat head with eight massage points. It brightens and revitalizes, stimulates collagen production, and uses heat to encourage cell regeneration – firming and strengthening skin structure.



Exclusive Series

Massage direction

When using Neno Beauty Algo, move the head from the center of the face outward and gently upward, following the guide lines in the illustration. For optimal results, comfort, and safety, always use with a conductive gel suitable for RF and EMS or a light, oil-free skincare product. Never use on dry skin – this may cause discomfort, irritation, and reduce effectiveness.



Exclusive Series

Memory function

No need to reset each time – Neno Beauty Algo remembers your last settings and starts with them automatically. One press is enough to enjoy a comfortable, personalized treatment.



Exclusive Series

Timer

Want to fully relax during skincare? The built-in timer guides the session for you: after three minutes, a gentle signal reminds you to switch sides of the face, and after six minutes, it indicates the end. With Neno Beauty Algo, every treatment is balanced and effortless.

6 min



Exclusive Series

Safety and hygiene guidelines

- **Alternate heads:** Use the diamond and flat heads on different days. While repeating the same mode isn't harmful, doing so several times a week may cause irritation.
- **Keep clean:** After each use, wipe the device with a dry or slightly damp cloth with mild soap. For disinfection, use a wipe soaked in isopropyl alcohol (up to 70%). Avoid harsh detergents, bleach, or immersion in water.
- **Do not use on implants:** Avoid use if you have electronic implants (e.g. pacemaker, insulin pump) or metal implants in the head/neck (dental, bridges, crowns, braces, hearing implants).
- **Avoid sensitive areas:** Do not use directly under the eyes, on eyelids, or on the thyroid area. Around eyes, nose, and mouth, always use the lowest intensity.
- **Moderation:** Best results come from consistent care 2–4 times a week (for 6 minutes), with breaks and varied modes.
- **Not for pregnancy or breastfeeding:** Use is strictly contraindicated during pregnancy and breastfeeding.
- **Before use:** Read all safety precautions in the user manual.



Exclusive Series

Accessories included

The elegant set includes protective covers for both massage heads, a storage pouch, and a microfiber cloth. Thanks to these, using Neno Beauty Algo is not only effective but also hygienic and safe.



Exclusive Series

neno[®] beauty

ALGO

SPECIFICATION

Model:	Neno Beauty Algo
Power supply:	5V 0.75A
Power:	12W
Battery:	1200 mAh
Timer:	6 min
Material:	ABS+PC+PP+zinc alloy
Dimensions:	54.5×45.3×190 mm
EAN code:	5902479677187
Box dimensions:	148×245×60 mm
Box weight:	0.605 kg
PN:	NEN-DOM-PI004
Set contents:	Massager, storage pouch, microfiber cloth, type-C cable, user manual
Contact:	Tel.: 731 594 301; mail: info@neno.pl
Producer:	KGK Trend sp. z o.o. Ujastek 5B, 31-752 Cracow



AWARDS FOR
KGK TREND
The Owner of
Neno brand



If you don't know which product to choose, need help, support, or have questions for our experts, call or write directly to the Neno Customer Support Center. Helpline: 731 594 301, mail: info@neno.pl