

neno[®] beauty

NAYA

FACE MASSAGER
WITH EMS TECHNOLOGY
AND THERMOTHERAPY FUNCTION

Exclusive Series



EMS technology

Exclusive Series

Discover the Neno Beauty series

The Neno Beauty series was created to make everyday skincare simple and enjoyable. In the comfort of your own home, you can enjoy the effects of intensive care – your skin regains energy, radiance, and a healthy look.



Exclusive Series

Discover Neno Beauty Naya

Neno Beauty Naya combines three modern technologies – EMS, thermotherapy, and vibration massage – which support the natural firming and brightening of the skin. Thanks to its dual massage head and five intensity levels, you can easily adjust the performance to your needs and create your own skincare ritual.

EMS

Electrical Muscle Stimulation: stimulates facial muscles, improving their tone and skin elasticity

Double massage head

Designed for precise stimulation of facial muscles

Thermotherapy

Gently warms the skin to 42°C, promoting microcirculation and increasing its elasticity.

Vibration massage

Reduces swelling and brings a feeling of relaxation

Timer

Automatic 6-minute timer and notification of massage area change every 15 seconds



Memory

Memory function

Easy to use

One function button

5 levels

5 intensity levels

Wireless

Wireless device, compact design

Stand

Ensures stable and hygienic storage of the massager

Exclusive Series



Thermotherapy



Vibration massage



EMS



Double head



Memory function



5 levels




Wireless



Timer



Exclusive Series




EMS ***technology***

EMS (Electrical Muscle Stimulation) uses gentle electrical impulses that stimulate muscle fibers to perform short, controlled contractions. This improves skin firmness and elasticity, making facial features appear more defined. During the treatment, you may feel a slight tingling sensation – this is a natural reaction of the skin to EMS impulses. Regular short sessions deliver the best results, with effectiveness depending on consistency and individual predispositions.



Exclusive Series




Thermotherapy function

Thermotherapy uses heat to improve circulation and support the skin's natural regenerative processes. Dilated blood vessels transport oxygen and nutrients more effectively, while muscles relax after a long day. This is not only skincare, but also a moment of pleasant relaxation for you.

Exclusive Series

Vibration massage

Vibration massage is based on quick, rhythmic pulsations that introduce gentle vibrations to the facial skin. This type of massage stimulates microcirculation, relieves tension, and gives the complexion a radiant appearance. It's a simple way to feel relaxed in just a few minutes.



Exclusive Series

Dual massage head

The ergonomic dual massage head was designed for precise stimulation of facial muscles. It naturally adapts to the contours of the face – from the jawline, through the cheeks, to the neck. This ensures the treatment is smooth, comfortable, and effective.




Exclusive Series

A close-up photograph of a woman's face. She is holding a white, handheld facial device with two brown rollers against her cheeks. She has a slight smile and is looking towards the camera. The background is a solid light blue.

5

intensity levels


With 5 intensity levels, you can adjust the device to your skin's sensitivity and desired effects. It is recommended to start at lower levels and gradually increase as you continue treatments – this helps the skin adapt and makes the experience more comfortable.

A dark brown horizontal bar at the bottom of the page. It features five bright, glowing white circles spaced evenly across its width.

Exclusive Series

Massage directions

When using Neno Beauty Naya, move the head from the center of the face outward and gently upward following the guide lines in the illustration. For optimal results, comfort, and safety, use the device only with a conductive gel suitable for RF and EMS technology, or with a light, oil-free facial cosmetic. Never use the device on dry skin – this may cause discomfort, irritation, and reduce treatment effectiveness.



prohibited area

recommended
areas of low usage



Exclusive Series



Safety and hygiene guidelines

- Maintain hygiene – after each use, clean the device by wiping it with a dry cloth or one slightly dampened with warm water and mild soap. For disinfection, wipe the heads with a cloth soaked in isopropyl alcohol solution (up to 70%). Avoid aggressive detergents, bleach, and immersing the device in water.
- Do not use if you have any electronic implants (e.g. pacemaker, insulin pump) or metal implants in the head and neck area (including dental implants, bridges, crowns, braces, or metal hearing implants).
- Avoid sensitive areas – do not use directly under the eyes, on the eyelids, or on the thyroid area of the neck. In the eye, nose, and mouth areas, always use the lowest intensity level.
- Regular but moderate use – best results come from consistent care 2–4 times per week for 6 minutes, with breaks and varied modes.
- The device must not be used during pregnancy or breastfeeding.
- Before using, read all precautions included in the user manual.

A close-up portrait of a woman with dark hair pulled back, holding a white, ribbed, handheld facial device with two silver spherical electrodes against her cheek. She has a soft smile and is looking directly at the camera. The background is a solid, light beige color.

Memory function

Neno Beauty Naya remembers the last settings used, so when you turn it on, it returns to your preferences without reconfiguration – quick, convenient, and consistent with your routine.

Exclusive Series

Timer

The timer in Neno Beauty Naya guides the rhythm of the treatment: 6 minutes total, with a signal every 15 seconds to remind you to change the area. After the time has elapsed, the device switches off automatically.



Exclusive Series


Accessories included

The transparent Neno Beauty Naya stand ensures stable storage of the product. You can place the device on your dressing table, bathroom shelf, or bedside table – always within reach, ready to use, and stored hygienically and elegantly.

The set also includes a soft microfiber cloth for daily cleaning and a USB cable for quick and convenient charging.



USB-C
cable



microfiber
cloth

stand



Exclusive Series

нено[®] beauty

NAYA

TECHNICAL DATA

Model:	Neno Beauty Naya
Power supply:	5V 0.7A
Power:	6W
Battery:	1000mAh
Timer:	6min
Material:	ABS + PC + silicone + zinc alloy
Dimensions:	84x45x68mm
Weight:	172g
EAN code:	5902479677163
Box dimensions:	147x207x60 mm
Box weight:	0,53kg
PN:	NEN-DOM-PI002
Set contents:	massager, stand, microfibre cloth, Type-C cable, user manual
Contact:	Tel.: 731 594 301; mail: info@neno.pl
Producer:	KGK Trend sp. z o.o. ul. Ujastek 5B, 31-752 Cracow



AWARDS FOR
KGK TREND COMPANY
Owner of the Neno brand:



If you do not know which product to choose, need help or support, or have questions for our experts, please call or write directly to the Neno Customer Support Centre.
Hotline: 731 594 301, mail: info@neno.pl