

bela ELECTRONIC BREAST PUMP



HYGIENE AND COMFORT



15 DIFFERENT MODES



FOR DELICATE SKIN















EASY STERILIZATION

The parts of the breast pump have to be regularly cleaned and sterilized. You can use boiling water or a steam sterilizer, also used to disinfect baby bottles.









UNIVERSAL

The Neno Bella breast pump can be used by any new mother, regardless of her breast size.









TECHNICAL SPECIFICATION

Model: neno Bella

Weight:	417g
Dimensions:	105x53 mm
Charging port:	USB-C
Battery type and capacity:	3.7v 1400 mAh Li-lon
Modes of operation:	6 stimulation levels, 9 expression levels
Charging time:	2 hours
Automatic power off:	after 30 min
Suction in expression mode:	18,9 - 42,9 kPa
Suction in stimulation mode:	3,9-38,9 kPa
Package contents:	Breast pump, USB-C cable, milk bottle, suction cup, suction bowl, bottle stand, bottle teat, valve x2, user's manual

