neno-

inspired by **Children**, designed for **Parents**.





COMFORT AND HYGIENE



SMART BREAST PUMP



HIGH-CLASS SILICONE





















EASY STERILIZATION

Parts of the breast pump has to be cleaned and sterilized regularly. You can do this using boling water or a steam sterilizer, like the one used to disinfect baby bottles.





TECHNICAL SPECIFICATION

Model: neno Due

| Weight: | 785g |
|----------------------------|---|
| Size: | 81x66x62 mm |
| Battery capacity and type: | 3.7v 1400 mAh Li-Ion |
| Charging port: | USB-C |
| Modes of operation: | 5 stimulation levels, 9 expression levels |
| Charging time: | 2 hours |
| Automatic power off: | after 30 min |
| Milk expression suction: | 20,3 - 42,9 kPa |
| Stimulation mode suction: | 14,5 - 22 kPa |
| Package contents: | Breast pump, 2x milk bottle, 2x silicone tube, 4x valve, 2x breast shield cover, USB-C cable, user's manual |

