

Medic TO2 INFRARED THERMOMETER









BY Xblitz

just one click

NENO MEDIC TO2 CAN BE USED TO SAFELY AND EASILY MEASURE A BABY'S BODY TEMPERATURE. JUST AIM THE DEVICE SENSOR AT YOUR CHILD'S FOREHEAD AND PRESS THE BUTTON TO GET THE MEASUREMENT DISPLAYED ON THE SCREEN.

11

Èn°.

infrared technology

NENO MEDIC TO2 USES INFRARED TECHNOLOGY. THE DEVICE RECEIVES INFRARED LIGHT WAVES AND USES THEM TO DETERMINE THE BODY TEMPERATURE WITHOUT TOUCHING THE PATIENT.



modern design

ONE BUTTON, EASY TO READ DISPLAY SCREEN AND LIGHT CASING MAKE THE NENO MEDIC TO2 THERMOMETER VERY EASY TO USE.

quick measurement

USING THE **NENO MEDIC T02** THERMOMETER IS INCREDIBLY SIMPLE. JUST PLACE THE DEVICE SENSOR 1,5-5 CM AWAY FROM THE PATIENT'S FOREHEAD AND PRESS THE BUTTON. **NENO MEDIC T02** WILL MAKE THE MEASUREMENT IN THE BLINK OF AN EYE.

nenomedic t02 INFRARED THERMOMETER

| Measurement spot: | Forehead |
|------------------------|---|
| Measurement distance: | 1,5 - 5cm |
| Measurement range: | Body temperature 32.0°C ~ 42.2°C (89.6°F ~ 108°F) |
| Measurement accuracy: | ± 0.2°C (35.0°C - 42.0°C)/ 0.4°F (95.0°F - 107.6°F) |
| Measurement units: | Degrees Celscius (°C) or degrees Fahrenheit (°F) |
| Memory: | 32 measurements |
| Temperature alert: | Yes, at 38°C, alarm noise |
| Power supply: | DC3V Baterie AAAx2 - not included |
| Automatic power-off: | After 1 min of inactivity |
| Operating temperature: | 16.0°C ~ 35.0°C |
| Relative humidity: | 15 - 80% |
| Storage temperature: | -20°C ~ 55°C |
| Storage humidity: | 15 - 93% |
| Dimensions: | 150x36x18 mm |
| Weight: | 51g |
| | |

