Meno* vapore

portable mesh nebulizer

♥ Every Mom

is a Superhero!







aerosol particle size













3 sizes of masks


3 speed modes

cordless











molecules MMAD

The effectiveness of inhalation depends on the particle size of the drug. Neno Vapore produces particles with an MMAD of $<5\mu m$ so that the drug reaches deep into the lower airways.












Large medication container

Thanks to the 6 ml medication container, nebulisation can be carried out without frequent refilling. This saves time and ensures continuity of therapy.















SPECIFICATION:

Power supply: DC 3.7 \

Battery: Lithium-ion 750 mAl

Power consumption: <4.0 V MMAD particle size: <5 μm

Nebulization speed: 0.15 ml/min - 0.90 ml/min (3 modes

Medicine tank capacity: 6 m

Continuous operation mode: Permanent or variable (cycle 10 minutes)

Respirable fraction: 65% Main Unit Weight: 72 g

Main unit dimensions: 7.25x6.75x5.05 cm
EAN code: 5902479675053
PN code: NEN-MAM-INH005

Box weight: 370 g

Box dimensions: 17x13.5x8.5 cm



AWARDS FOR KGK TREND COMPANY

The owner of Nho brand









