

neno-Freedom Geries

The Freedom Series breast pumps are designed for women who value convenience and independence above all else. Lightweight and portable, they change the face of breast pumping by being hands-free. This allows any mother to use her hands for other activities while expressing breast milk for her baby.





Pearl in Neno's crown

The elegant Perla is the first such breast pump in the Neno range. It combines an innovative hands-free formula with the functionality of a classic bottle pump.

This brings a whole new dimension to everyday breastfeeding.







Stimulation phase

Fast, intense pulsation. Stimulates milk flow. The breast pump reproduces the baby's natural sucking action.

Deep extraction phase

Gentle pumping to express milk from the breast.

Mixed phase

An innovative phase that has a direct effect on stimulating lactation and helps with complex breast emptying and lack of milk.

Three modes of operation

You can use the breast pump in three modes. In each of them, you can adjust the suction power in a range of 9 levels.













Practical container


Don't worry if you need to express your breastmilk but don't have a bottle or pouch with you. The breast pump container with lid can also be used for storage.





neno Premium Simple function panel

The Neno Perla breast pump touchpad provides comfortable and intuitive operation. You can easily select modes and adjust the suction power, matching it perfectly to your current needs.















Comfortable use

After expressing the milk, you can easily store it in the fridge by securing the container with the lid. You can also pour it straight into a bottle!



Durable battery

The Neno Perla is equipped with a high-capacity battery, eliminating the need for frequent recharging. This is a great convenience for women with active lifestyles who want to pump regularly at the same time.





quick sterilisation Maintaining good hygiene when

Maintaining good hygiene when expressing milk reduces the risk of your baby coming into contact with harmful organisms. Pay particular attention to frequent sterilisation of parts that come into direct contact with milk.









choose

FOR ONE OR TWO BREASTS











SPECIFICATION

Model: Neno Perla

Work modes: Stimulation phase L1-L9

Deep expression phase L1-L9

Mixed phase L1-L9

Capacity of milk container: 120 ml

Acumulator: Li-ion 3.7 V 1100 mAh

Charging time: about 2 hours

Operating time: up to 120 minutes

Input: 5V 1A

Charging cable: USB-C

Dimensions: 13x11x6,5 cm

Weight: 277 *g*

EAN code: 5902479675565

Package dimensions: 175 x 165 x 80 mm

Package weight: 826 g

PN: NEN-MAM-LK015

Included: Pumping station • USB-C cable • Silicone breast shield 26 mm

• 2× Valve • Milk container • Milk container cap

• Breast pump connector • Velvet pouch • Silicone brush

Manual

neno Premium

Awards for KGK TREND company
Owner of Neno brand:







