

neno[®] atemo

Medical Device

PORTABLE
ULTRASONIC
NEBULIZER

♥ *Every Mom
is a Superhero!*



2 ways of
nebulization

MMAD
<5 μ m

fine
particles



for children
and adults



cordless



2 nebulization modes

The nebulizer offers two modes of medication delivery – through a mask or a mouthpiece. This allows you to choose the inhalation method that best suits your preferences and situation. Such flexibility makes therapy convenient and perfectly adapted to your daily routine.



Medical device

Our nebulizer is a certified medical device that meets the required quality standards. By choosing Neno Atemo, you're opting for a reliable product that ensures both safety for your family and effective performance.





Wireless operation

No more tangled cables or searching for a power outlet. Neno Atemo is a fully wireless nebulizer that guarantees comfort and convenience wherever you are. Whether you plan to do your inhalation at home, during a business trip, or elsewhere — you can use it anytime, anywhere. Feel the difference with Neno Atemo.

Fast inhalation

Neno Atemo is ready to use in seconds. Its efficient system quickly delivers medication to the respiratory tract, bringing relief and reducing unpleasant symptoms. It's the perfect solution both during cold and flu season and for everyday respiratory care.



For children and adults

The Neno Atemo nebulizer is designed for both adults and children, making it suitable for the whole family. It can be used not only during illness but also for regular inhalation therapy as part of your daily routine.



Compact size

Neno Atemo stands out for its compact design, making it easy to take on trips, to work, or even for a walk. Despite its small size, it delivers powerful and effective therapy wherever you are.



Quiet operation

Our nebulizer operates extremely quietly, allowing for use anytime — day or night. It won't disturb others, making inhalation calm, comfortable, and stress-free.





Perfect for travel

Thanks to its small and lightweight design, Neno Atemo is the ideal travel companion. It fits easily into your luggage, making it perfect for family trips and business travel alike.

Automatic shut-off

The nebulizer features an automatic shut-off function — it stops working after 20 minutes of continuous operation. This ensures safe and convenient use. The intermittent mode (20 minutes on, 10 minutes off) allows you to tailor therapy time to your needs.

A young girl with brown hair is sitting on a beige couch, wearing a white nebulizer mask with blue straps. She is holding an open book. A teddy bear is visible behind her. In the foreground, there is a circular graphic with a teal arrow pointing clockwise and the text '20 MIN' in teal.

**20
MIN**

The image features a white Neno Atemo nebulizer on the left, with its clear plastic chamber open, revealing a fine mist of particles. To the right, a detailed molecular model is shown, consisting of dark spheres connected by thin rods, representing a complex chemical structure. The background is a soft, light blue gradient.

particles
MMAD **<5**
 μm

Fine particles are essential for effective therapy. Neno Atemo produces aerosol particles with a mean diameter of less than 5 μm (MMAD <5 μm), allowing the medication to reach deep into the lower respiratory tract for fast and efficient results.

Long battery life

Equipped with a powerful and durable lithium battery, Neno Atemo provides up to 2.5 hours of continuous use on a single charge. You can perform inhalations freely without worrying about unexpected power loss.



Large medication container

The large 8 ml medication container ensures comfortable and uninterrupted inhalation sessions, reducing the need for refills. Enjoy longer and smoother therapy without unnecessary pauses.



Complet set

In addition to the nebulizer itself, the set includes two masks (for children and adults), a mouthpiece, a USB-C cable, a storage pouch, and a user manual. Everything you need for a complete and comfortable therapy experience with Neno Atemo.



SPECIFICATION

Power supply:	Lithium battery DC 3.7V or external power supply 5V 1A (via USB-C)
Power consumption:	<2VA
MMAD particle size:	<5µm
Nebulization rate:	0.2ml/min
Medication cup capacity:	8ml
Ultrasonic frequency:	110 ±10kHz
Medication residue:	<1ml
Operating environment:	10°C–40°C (50–104°F), 10%RH–95%RH, 86kPa–106kPa
Automatic shut-off:	Yes, after 20 minutes
Noise level:	≤50dB
Operating time:	Up to 2.5h
Charging time:	Up to 1.5h
Weight:	46g
Dimensions of the main unit with medication cup and mesh:	76 × 56 × 47mm
Box weight:	167g
Box dimensions:	107 × 125 × 57mm
EAN code:	5902479677750
PN:	NEN-MAM-INH007
Set contents:	Nebulizer with medication cup, adult mask, child mask, mouthpiece, storage pouch, USB-C charging cable, user manual
Contact:	+48 731 594 301, E-mail: info@neno.pl
Importer:	KGK Trend sp. z o.o. Ujastek 5B, 31-752 Cracow, Poland

If you're not sure which product to choose, need help or support, or have questions for our Neno midwives or lactation experts, call or write directly to the Neno Customer Support Center.
Helpline: +48 731 594 301, Email: info@neno.pl

AWARDS FOR
KGK TREND
COMPANY
Owner of the
Neno brand:

