

Touchless thermometers

Model	 Neno Medic T05	 Neno Medic T02	 Neno Medic T07	 Neno Medic T06
Measurement distance	1.5–5cm	1.5–5cm	1–3cm	1–5 cm
Measurement accuracy	±0.2°C in range 35~42°C, (±0.3°C in other ranges)	± 0.2°C in range 35°~42°C (± 0.3°C in other ranges)	forehead ±0.2°C in range 36~39°C, (±0.3°C in other ranges) ear ±0.2°C in range 36~39°C, (±0.3°C in other ranges) object/room: ±1°C	forehead ±0.2°C in range 36~39°C, ±0.3°C in ranges 22~36°C and 39~43°C object ±1.0°C/±2.0°F
Measurement range	34~43°C (forehead) 0~100°C (object)	32.0°~42.2°C (forehead)	22~43°C (forehead) 34~43°C (ear) 0~100°C (object)	22~43°C (forehead) 0~100°C (object)
Measurement place	Forehead or object	Forehead	Forehead, ear or object	Forehead or object
Units of measurement	Stopnie Celsjusza (°C) lub Fahrenheita (°F)	Degrees Celsius (°C) or Fahrenheit (°F)	Degrees Celsius (°C) or Fahrenheit (°F)	Degrees Celsius (°C) or Fahrenheit (°F)
Measurement memory	34	32	20	20
Temperature alert	✓	✓	✓	✓
Power Supply (Batteries)	2×AAA	2×AAA	2×AAA	2×AAA
Automatic switch-off	after 10 seconds of inactivity	after 60 seconds of inactivity	after 10 seconds of inactivity	after 10 seconds of inactivity
Additional bag included	×	×	✓	✓
Dimensions	150×37×16.3 mm	150×37×16.3 mm	145.5×45.4×43.4 mm	149.8×47.3×39.1 mm
Weight (without batteries)	105 g	51 g	72 g	73 g
Included	Thermometer, User manual	Thermometer, User manual	Thermometer, Bag, 2×Batteries AAA, User manual	Thermometer, Bag, 2×Batteries AAA, User manual