

# neno<sup>®</sup> pressione

BLOOD PRESSURE MONITOR



Polish  
voice  
commands



Arrhythmia  
detection

22-42cm



Wide cuff  
range



Medical  
device

AWARDS FOR  
KGK TREND COMPANY  
The owner of Neno brand:



# arrhythmia *detection*



The Neno Pressione has an intelligent system that generates an alert when an irregular heartbeat occurs. The detection of an anomaly is indicated by the ECG curve symbol, located in the lower left corner of the display. This option allows you to be aware in real time that a doctor's appointment should be made.

**Please note that blood pressure monitors are intended for self-monitoring and not for diagnosis. If you experience any disturbing symptoms, consult your doctor immediately.**

# wide *cuff range*

The electronic blood pressure monitor ensures comfortable use thanks to its cuff with a wide adjustment range (22 cm - 42 cm). The flexible cuff, made of high-quality materials, allows you to adapt the measurement to your individual needs.




**22 - 42 cm**

clear  
*LCD display*



The clear and easy-to-read display serves for the correct and safe use of the blood pressure monitor. Through it, you can access the numerical result of the measurement and the scale of correctness of pressure expressed in colours.

Use  
*wherever you like*



The Neno Pressione is equipped with two types of power supply. Traditional mains power via a USB C socket works well at home. Battery power enables you to use the blood pressure monitor wherever you go.



# comfortable *cover*



*Quick  
zip*



*Practical  
handle*



The included case is the best way to store your blood pressure monitor. The Neno Pressione will always have a place not only in your heart, but also in the case.



The included carrying case makes it easy to safely transport and store the device when travelling, at doctor's appointments, hospital stays and also on a day-to-day basis.

# *TRAVEL* companion



# simple operation



The operation of the blood pressure monitor is so intuitive that both younger and older users can easily handle it. Measurements are taken at the press of a button. In addition, a quick reference guide containing the most important information is included.





# *memorising results!*



# 2X 99 stored results

The result memory function is particularly useful for people who need to monitor their heart rate on an ongoing basis. The Neno Pressione upper arm blood pressure monitor memorises 99 results for each of two people. You can read out the memorised data at any time, allowing you to track differences between measurements.



mmHg



kPa



choice *of unit*

Depending on your choice, the results are expressed in either kilopascals or millimetres of the column of mercury.

# unique design



The device attracts attention for its build quality and unique design. The glossy coating in an elegant black colour makes the Neno Pressione stand out from other electronic blood pressure monitors on the market.


# The result at hand

The Neno Pressione is a device that will interpret your result. Next to the pressure data, glowing colours will appear to suggest how you should read the measurement.

result within normal range



result out of the norm



incorrect result



# Medical *product*



The Neno Pressione is an electronic upper arm blood pressure monitor classified as a medical device, which is why you will find it in pharmacies, specialist shops and partner shops.



# nenopressione

BLOOD PRESSURE MONITOR

<b>LED display:</b>	98×83 mm
<b>Measurement method:</b>	oscillometric
<b>Measurement range:</b>	for a pressure of 0-295 mm Hg (0-39.3 kPa) for a pulse of 40 to 199 beats per minute
<b>Units of measurement:</b>	millimetres of mercury (mmHg) or kilopascals (kPa)
<b>Power supply:</b>	battery operated (4× AA battery) or mains (USB-C DC 5V socket)
<b>Memory:</b>	up to 99 measurements for 2 users
<b>Cuff:</b>	Shoulder 22-42 cm
<b>Working conditions:</b>	temperature: 5-50°C air humidity 15-85% atmospheric pressure 70-106 kPa
<b>Storage conditions:</b>	temperature -20-55°C air humidity 10-93% atmospheric pressure 70-106 kPa
<b>Dimensions:</b>	141×120×54 mm
<b>Weight:</b>	304 g
<b>EAN:</b>	5902479673936
<b>CN:</b>	90262020
<b>PN:</b>	NEN-MAM-CI001



If you don't know which product to choose, need help, support or have questions for our midwives, neno lactation experts, call or email the neno Customer Support Centre directly. Helpline: 731 594 301 mail:info@nenopressione.pl

AWARDS FOR  
KGG TREND COMPANY  
The owner of Neno brand:



MADE IN POLAND  
2022